



## MediMix

Combined Mixing and Medicating System

# MediMix

## MediMix™ Combined Mixing and Medicating System

Many types of additives, medicines, and supplements are being used in livestock applications.

All of the above, most of the time, are injected into the drinking water system by a hydraulic proportioning medicators.

Part of the additives attend to separate and to sink down inside the additive container.

It is highly recommended to allow proper agitation of the additives during the entire medicating process, in order to assure fully homogenized solution.

TEFEN is offering to the growers a unique, high end system of medication and agitation.

The MediMix™ is offering:

1. Quality proportioning hydraulic medicator.
2. Electric overload mechanism for long lasting
  - 2.1. Main switch
  - 2.2. External selector switch 3 modes:  
Manual / closed / automatic (timer mode)
  - 2.3. On/Off timer to control working hours/pause
  - 2.4. Contactor (relay)
  - 2.5. Water protected three-phase plug
3. Rigid quality 200 liter PE tank.
4. Galvanized metal frame.
5. State of the art SS agitator with an advanced motor and gear.

\* Dimensions: length 652 mm, width 652 mm, hight 882 mm.

\* weight: approximately 27 kg.

